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The Fourier Transform InfraRed spectrophotometer measures an infrared spectrum by Fourier-transform of an interferogram.

The Physics of Fourier Transform Spectrophotometry

An FTIR uses one of several optical systems; the IRAffinity-1 relies on a Michelson interferometer (Fig. 1.1). After passing through the aperture, light is turned into a parallel beam by the collimator mirror and enters the beam splitter. A germanium film, deposited on a potassium bromide substrate via evaporation, comprises the beam splitter; it splits the single beam into two, reflecting one to the fixed mirror and transmitting the other to the moving mirror. Both mirrors reflect their beams back to the beam splitter; part of each returning beam is reflected and transmitted. The transmitted light from the fixed mirror and the reflected light from the moving mirror recombine and interfere with each other as they travel towards the collecting mirror. The interference is either constructive or destructive.


Fig. 1.1 Michelson Interferometer

Assume that the light source emits monochromatic light of wavelength λ (cm). When the distance l_1 between the fixed mirror and the beam splitter is equal to the distance l_2 between the moving mirror and the beam splitter, the optical path difference between the two beams, $\chi = 2(l_1 - l_2)$, is equal to zero, and the beams are in phase. While in phase, the beams interfere constructively with each other (Fig. 1.2 A, B). As the moving mirror is displaced $\lambda/4$ cm, the optical path difference becomes $\lambda/2$ cm, and the two beams are out of phase, interfering destructively (Fig. 1.2 A, C). Thus, the two beams interfere constructively with each other when $\chi = n\lambda$ and destructively when $\chi = (n + 1/2)\lambda$ where n is an integer.


Fig. 1.2 Interference

Equation 1.1, extrapolated from the above principles, calculates the intensity $I^*(\chi)$ of light (wavelength λ) incident to the detector.

$$I^*(\chi) = 4RTS(\lambda) \left[\frac{1}{2} + \frac{1}{2} \cos 2\pi \frac{\chi}{\lambda} \right] \dots\dots\dots (1.1)$$

- where R: energy reflected by the beam splitter
- T: energy transmitted by the beam splitter
- S (λ): radiation energy from the light source

The intensity of the light observed by the detector is a function of Equation 1.1. $I(\chi)$ denotes the light intensity, and the wavenumber σ (cm^{-1}) replaces the wavelength λ .

$$I(\chi) = 4RTS(\lambda) \cdot \frac{1}{2} \cos 2\pi \sigma \chi$$

$$= B(\sigma) \cos 2\pi \sigma \chi \dots\dots\dots (1.2)$$

where $B(\sigma) = 4RTS(\lambda) \cdot 1/2$

The signal $I(\chi)$ observed by the detector is called an interferogram, and $4RT$ is labeled beam splitter efficiency. If polychromatic light is emitted instead of monochromatic light, $I(\chi)$ is given by the integration of 1.2 with respect to wavenumber.

$$I(\chi) = \int_0^\infty B(\sigma) \cos 2\pi \sigma \chi \delta \sigma \dots\dots\dots (1.3)$$

Equation 1.3 demonstrates that $I(\chi)$ is a Fourier cosine transform of spectrum $B(\sigma)$. Thus, an inverse Fourier cosine transform of $I(\chi)$ recovers the original spectrum $B(\sigma)$.

$$B(\sigma) = \int_{-\infty}^{+\infty} I(\chi) \cos 2\pi\sigma\chi \delta\chi \dots\dots\dots (1.4)$$

While the conventional dispersive spectrometer directly determines the intensity of $B(\sigma)$ at certain wavenumbers, the FTIR detector observes the interferogram $I(\chi)$ which must be Fourier transformed to obtain the spectrum.

Although Equation 1.4 maintains that the interferogram should be recorded from $-\infty$ to $+\infty$, such limits are impractical. In reality, the integration is limited to an optical path difference L , determined by the range of the moving mirror. The experimental spectrum $B'(\sigma)$ differs from the theoretical spectrum $B(\sigma)$ accordingly.

$$B'(\sigma) = \int_{-L}^L I(\chi) \cos 2\pi\sigma\chi d\chi$$

$$= \int_{-\infty}^{+\infty} A(\chi) I(\chi) \cos 2\pi\sigma\chi d\chi \dots\dots\dots (1.5)$$

where $A(\chi) = \begin{cases} 1 & \text{if } |\chi| \leq L \\ 0 & \text{if } |\chi| > L \end{cases} \dots\dots\dots (1.6)$

According to Convolution Theory,

$$B'(\sigma) = B(\sigma) * F(\sigma) \dots\dots\dots (1.7)$$

Where $F(\sigma)$ is a Fourier transform of the function $A(\chi)$ in Equation 1.6.

$$F(\sigma) = 2L \sin(2\pi\sigma L) / (2\pi\sigma L)$$

$$= 2L \text{sinc}(2\sigma L) \dots\dots\dots (1.8)$$

The measured spectrum $B'(\sigma)$ is a convolution of theoretical spectrum $B(\sigma)$ with $F(\sigma)$. $F(\sigma)$, a Fourier transform of $A(\chi)$, is called instrument function or instrumental line shape ILS.

Functions $A(\chi)$ and $F(\sigma)$, represented by Equations 1.6 and 1.8 respectively, are given in Fig. 1.3 (a). The function $F(\sigma)$, as defined by Equation 1.8, first becomes zero at $(= \pm (1/2L))$, where the half-width value is $0.605/L$. The greater the optical path difference L , the smaller the half-value width; when the moving mirror is driven further, the resolution improves.

The instrument function $F(\sigma)$ involves marked submaxima, called side-lobes [Fig. 1.3 (a)]. If $A(\chi)$, as defined by 1.6, is replaced by the triangular function defined by 1.9, Equation 1.10 gives its Fourier transform.

$$A(\chi) = \begin{cases} 1 - |\chi|/L & \text{for } |\chi| \leq L \\ 0 & \text{for } |\chi| > L \end{cases} \dots\dots\dots (1.9)$$

$$F(\sigma) = L \sin^2(\pi\sigma L) / (\pi\sigma L)^2 = L \text{sinc}^2(\sigma L) \dots\dots\dots (1.10)$$

Fig. 1.3 (b) illustrates that the instrument function defined by 1.10 has worse separation, but smaller submaxima, than that defined by 1.8. Using $A(\chi)$ to reduce the instrument function submaxima called apodization; $A(\chi)$ is called the apodization function.


Fig. 1.3 Apodization with "Rectangular" and "Triangular" Functions

Table 1.1 lists the IRAffinity-1 apodization parameters and their characteristics.

Table 1.1 Apodization Function

Parameter	Apodization Function	Instrument Function	Half-width value	Submaxima residue
None, Box Car	$A(\chi) = 1$	$2L \operatorname{sinc}(2\pi\sigma L)$	0.605/L	-21%
Triang	$A(\chi) = 1 - \frac{ \chi }{L}$	$L \operatorname{sinc}^2(\pi\sigma L)$	0.88/L	+4.5%
Sqr. Triang	$A(\chi) = \left(1 - \frac{ \chi }{L}\right)^2$	$\left\{\frac{4L}{(2\pi\sigma L)^2}\right\} 1 - \operatorname{sinc}(2\pi\sigma L)$	1.18/L	0.7%
Bessel	$A(\chi) = \left\{1 - \left(\frac{\chi}{L}\right)^2\right\}^2$	$L(2\pi \times \sigma)^{5/2} J_{5/2}(2\pi \times \sigma)$	0.952/L	-4.1%
Cos	$A(\chi) = \frac{1}{2} \left\{1 + \cos\left(\frac{\pi \cdot \chi}{L}\right)\right\}$	$\frac{\operatorname{sinc}(2\pi\sigma L)}{2\pi\sigma(1 - 4L^2\sigma^2)}$	1.00/L	-2.7%
Happ	$A(\chi) = 0.54 + 0.46 \cos\left(\frac{\pi\chi}{L}\right)$	$\left\{\frac{0.54}{\pi\sigma} + \frac{(0.46) \cdot 4\pi\sigma L^2}{\pi^2 - (2\pi\sigma L)^2}\right\} \operatorname{sinc}(2\pi\sigma L)$	0.91/L	-0.6%

Notes: (a) $A(\chi) = 0$ ($|\chi| > L$)

(b) Ratio of the maximum submaxima peak size to the center peak size (shown in %).

The information in Table 1.1 supports using the "None" or "Box Car" ("Rectangular" Function) parameter for samples requiring high resolution, such as gas, and the "Happ" ("Happ-Genzel" Function) or "Sqr. Triang" ("Sqr-Triang" Function) parameter, which has fewer submaxima peaks, when high resolution is unnecessary.